

## VELP Scientifica Automatic Digestion Units



**DKL 8**



**DKL 12**



**DKL 20**



**DKL 42/26**

Number of positions		8 x 250 ml	12 x 250/400 ml	20 x 250 ml	42 x 100 ml
Lift movement		Automatic	Automatic	Automatic	Automatic
Temperature range		Ambient to 450 °C	Ambient to 450 °C	Ambient to 450 °C	Ambient to 450 °C
Stability and precision of heating block temperature		± 0.5 °C	± 0.5 °C	± 0.5 °C	± 0.5 °C
Program library		30 standard + 24 user-programmable	30 standard + 24 user-programmable	30 standard + 24 user-programmable	30 standard + 24 user-programmable
Selectable ramps		up to 4 per program	up to 4 per program	up to 4 per program	up to 4 per program
Digestion time range		from 1 to 999 minutes	from 1 to 999 minutes	from 1 to 999 minutes	from 1 to 999 minutes
Supplied with	Lift	•	•	•	•
	Suction cap & drip tray	•	•	•	•
	Sample rack	•	•	•	•
	Test tubes	•	•	•	•
Dimensions (WxHxD) mm		210x690x540	266x690x540	322x690x584	322x690x584
	(in)	(8.3x27.2x21.3)	(10.5x27.2x21.3)	(12.7x27.2x23.0)	(12.7x27.2x23.0)
Weight	Kg	19.7	23.3	30.8	33.5
	(lb)	(43.5)	(51.4)	(68.0)	(74.0)
Power Supply		115 or 230 V	115 or 230 V	230 V	230 V
Power		1150 W	1500 W	2300 W	2300 W



[www.velp.com](http://www.velp.com)